

Actichlor™ Plus



Convenient and Effective Chlorine Tablets for Broad-spectrum Cleaning and Disinfection

Broad-spectrum cleaning and disinfection tablet for floors and general surfaces.

Actichlor Plus is a chlorine disinfectant product for all aspects of surface and environmental disinfection. Actichlor Plus offers:

- Effective chlorine disinfection in tablet format
- Broad-spectrum efficacy suitable for a wide variety of applications, including blood spills and sporicidal disinfection
- Easy-to-use format for single-step cleaning and disinfection

The Actichlor Plus tablet combines a chlorine compatible detergent with NaDCC in a single tablet format

- Full spectrum of efficacy
- Excellent cleaning and disinfection performance in one easy step



Actichlor™ Plus

Broad-spectrum cleaning and disinfection tablet for floors and general surfaces.

Product characteristics

- White tablet releasing active chlorine
- 1 tablet (3.22g) releases 1g active chlorine
- Hygroscopic (keep away from humidity)
- pH after dilution (from 200ppm - 10,000ppm/or 1,000 and 10,000) with tap water: 6 -7
- Clear colourless solution after dilution
- Working solution can be used up to 24 hours

Microbiological properties

Claim/Conditions	Concentration	Contact time
Dirty Conditions		
Bactericidal according to EN 13727/EN 16615	1000 ppm (1 tab/l)	5 min
Yeasticidal according to EN 13624/EN 16615	1000 ppm (1 tab/l)	5 min
Fungicidal according to EN 13624	10000 ppm (10 tab/l) 5000 ppm (5 tab/l)	5 min 15 min
Tuberculocidal/mycobactericidal according to EN 14348	5000 ppm (5 tab/l)	15 min
Virucidal according to EN 14476	1000 ppm (1 tab/l)	5 min
Sporicidal against C. difficile according to EN 17126	10000 ppm (10 tab/l) 5000 ppm (5 tab/l)	5 min 15 min
Sporicidal according to EN 17126	10000 ppm (10 tab/l)	15 min
Clean Conditions		
Bactericidal according to EN 13727 (200ppm, 1 min) / EN 16615 (200ppm, 10 min)	200 ppm (1 tab/5L)	10 min
Yeasticidal according to EN 13624 (100ppm; 1 min) / EN 16615 (200ppm, 10 min)	200 ppm (1 tab/5L)	10 min
Fungicidal according to EN 13624	5000 ppm (5 tab/l) 2000 ppm (2 tab/l)	5 min 15 min
Tuberculocidal according to EN 14348	200 ppm (1 tab/5L)	15 min
Mycobactericidal according to EN 14348	1000 ppm (1 tab/l)	15 min
Virucidal according to EN 14476	500 ppm (1 tab/2l)	5 min
Sporicidal against C. difficile according to EN 17126	1000 ppm (1 tab/l)	5 min
Sporicidal according to EN 17126	10000 ppm (10 tab/l) 5000 ppm (5 tab/l)	5 min 10 min

Packaging and ordering information

Product	Packaging unit	Order code
Actichlor™ Plus	6 X 150 Tablets	3077170
Actichlor™ Dilution Bottle	1	10061468

ECOLAB EUROPE GmbH
Hofwiesenstrasse 349
CH-8050 Zürich
Switzerland

www.ecolab.com



ECOLAB DEUTSCHLAND GmbH
Ecolab-Allee 1
40789 Monheim am Rhein
Germany



Indications

Broad-spectrum cleaning and disinfection tablet for floors and general surfaces.

Instructions for use

Concentrate. Completely dissolve number of tablets according to required concentration in water (max. 30° C).

To dilute e. g. to a concentration of 1000ppm, dissolve one tablet in 1 litre of water; or for a solution of 10,000ppm, ten tablets in 1 litre of water.

Apply in sufficient quantity (+/- 30 ml/m²). Respect the indicated contact time and conditions for the required antimicrobial activity.

If used under clean conditions, clean surface upfront (use clean, fresh cloth/mop for each step).

Rinsing is not necessary unless the surface comes in contact with skin.

For the frequency of use and cleaning of the application equipment, refer to the hygiene plan in place.

Ingredients

100g contains 52,31g Troclosene sodium (CAS N° 2893-78-9). 1 tablet (3.22g) releases 1.0g available Chlorine.

Storage

Between 0°C and +25°C.

Biocide: use safely

Please read label and product information.



BIOCIDAL PRODUCT INTENDED FOR SURFACE
DISINFECTION (GROUP 1 - PT2)
FOR PROFESSIONAL USE ONLY

ECOLAB[®]
PROTECTING WHAT'S VITAL™