

FOR **EFFECTIVE NEUTRALISING** AFTER **ALKALINE CLEANING**

Sekumatic[™] FNZ

LIQUID NEUTRALISING AGENT

Based on citric acid, Sekumatic FNZ is the neutralising agent for use after alkaline cleaning.

Suitable for use with many alkaline cleaners, the tenside free product can be used in both thermal and chemo-thermal processes.

Sekumatic FNZ also provides:

- High biodegradability
- Ease of use
- Excellent material compatibility





UK Responsible Person Ecolab Ltd. Winnington Avenue, Northwich, Cheshire CW8 4DX United Kingdom



Sekumatic[™] FNZ

Liquid neutralising agent



5 L Canister

20 L Canister

WHY USE SEKUMATIC FNZ?

Suitable for use with high alkaline cleaners or when a hard water quality is used during the cleaning step, Sekumatic FNZ is an effective liquid neutraliser.

CONVENIENT FORMATS

Containing aqueous citric acid, Sekumatic FNZ is supplied in easy to handle 5L canisters or a choice of 20L canisters and a 200L drum.

The product can be conveniently integrated into existing procedures and protocols.

YOUR ADVANTAGES

Highly effective at a low dosage. Good material compatibility with acid resistant materials in medical areas.

CHEMO-PHYSICAL DATA

•	Density	20°C	1.16 g/ml
•	pH value* (concentrate)	20°C	Approx. 2.3
•	pH value* (solution, 1 ml/L in VE water)	20°C	2.8 - 3.2
•	Conductivity (solution, 1 - 3 ml/L in VE water)	20°C	2.2 - 4.6 mS/cm

*pH value measurements under operating conditions always using pH meter with temperature compensation.

PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Sekumatic FNZ Canister	2 x 5 L can	3023510
Sekumatic FNZ Canister	1 x 20 L can	3023500

INSTRUMENT REPROCESSING

APPLICATION AREA

For neutralisation after alkaline cleaning of laboratory glassware, surgical instruments, instruments composed of plastic or elastomers (e.g.: rubber) from anaesthesia, surgery and urology.

INSTRUCTIONS FOR USE

Program the machine according to manufacturer's instructions. Dose approx 1 ml/L using the machine's integral acid resistant dosing pumps.

INGREDIENTS

Aqueous citric acid solution.

STORAGE

Dangerous - respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: between +0°C and +25°C.

CERTIFICATES

On demand we are pleased to provide you with certificates.

MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.

ECOLAB HEALTHCARE Winnington Avenue, Northwich. Cheshire CW8 4DX United Kingdom



ECOLAB DEUTSCHLAND GmbH Ecolab-Allee 1 40789 Monheim am Rhein

Germany

