



## BROAD SPECTRUM READY-TO USE SURFACE CLEANING AND DISINFECTION WIPE

### Incidin™ OxyWipe S

#### BROAD SPECTRUM CLEANING AND DISINFECTION WIPE

The advanced formulation of Incidin™ OxyWipe S provides the benefits of a patented\* hydrogen peroxide solution with broad spectrum efficacy. It is effective against a wide range of pathogens and leaves no active residues.

- Active ingredient decomposing into water and oxygen
- Virucidal offering
- Excellent material compatibility
- Incidin OxyWipe S efficacy confirmed for liquid released from the wipe



\*European Patent EP1987121

# Incidin™ OxyWipe S

## BROAD SPECTRUM CLEANING AND DISINFECTION WIPE

### EFFECTIVE

- Virucidal disinfectant

### SAFE

- Active ingredient decomposes into water and oxygen
- Excellent material compatibility

### CONVENIENT

- Cleans and disinfects
- Ready to use formats

### MICROBIOLOGICAL PROPERTIES

EFFICACY	MEDICAL DEVICE	BIOCIDE (PT02)
Bactericidal according to EN 13727 (2min) / EN 16615 (2 min*)	2 min*	15 min*
Yeasticidal according to EN 13624 (15min) / EN 16615 (1min*)	15 min*	
Fungicidal according to EN 13624 (30min) / EN 16615 (method) (5min*)	30 min*	30 min*
Tuberculocidal according to EN 14348 (30min) / EN 16615 (method) (5min*)	30 min*	30 min* <sup>1</sup>
Mycobactericidal according to EN 14348 (30min) / EN 16615 (method) (5min*)	30 min*	
Effective against enveloped viruses according to EN 14476	30 sec*	-- <sup>2</sup>
Limited spectrum virucidal according to EN 14476 (Adenovirus - 30sec, Norovirus - 15min)	15 min*/**	-- <sup>2</sup>
Virucidal according to EN 14476	30 min**	30 min** <sup>2</sup>

\* dirty conditions \*\* clean conditions

<sup>1</sup> Mycobactericidal efficacy includes tuberculocidal efficacy.  
<sup>2</sup> Only full virucidal efficacy was applicable under EN 14476 at the time of dossier compilation and submission.

### MEDICAL DEVICE AND BIOCIDES USE SAFELY

Please read label and product information.

### INDICATIONS



**Medical Device:** Broad spectrum ready to use cleaning and disinfection wipe for non-invasive, invasive and non immersible medical devices.

**Biocide:** Disinfection of non-food contact surfaces in health care applications by wiping using impregnated RTU wipes.

### INSTRUCTIONS FOR USE

**Medical Device:**  
Take single wipe. Wet surface thoroughly, ensuring contact with entire surface. Let it dry and do not rinse. Do not reuse wipe: Disinfection can no longer be guaranteed. Close lid after use. Respect the indicated contact time and conditions for the required antimicrobial activity. If used under clean conditions, clean surface upfront (use new wipe for each step). The use of a cleaning and disinfection wipe does not replace regular cleaning. Do not use on surfaces sensitive to oxidative agents such as marble, copper or brass. For reprocessing of medical devices always observe manufacturers' instructions for use, incl. material compatibility. Terminal disinfection of semi-critical medical devices (i.e. probes & TEE probes): If immersion is not possible, clean the device and then disinfect by wiping as described above. Rinse thoroughly with water of drinking quality in this case after contact time elapsed.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

**Biocide:**  
Product is intended for cleaning and disinfection in one step. Always read the label or leaflet before use and follow all the instructions. When used under clean conditions: clean surface up front. Product should be applied to a dry surface. Wet surfaces completely. After use of the product allow surface to air dry. Do not rinse after use. Close container when not in use. Do not use wipes after dehydration. Throw away the container when empty. Do not use on surfaces sensitive to oxidative agents such as marble, copper or brass. Used wipes must be disposed in a closed container.

**Caution:** The product is classified as Biocide as well as Medical Device. Users must respect contact times and conditions for the intended use.

### COMPOSITION

100 g of the impregnation solution contains: 1.5 g Hydrogen peroxide (CAS N° 7722-84-1)

### STORAGE

Product should be stored in the original container between 5°C and 25°C in a cool, dry location out of direct sunlight. Period after opening: 3 months.



FOR PROFESSIONAL USE ONLY

# Incidin™ OxyWipe S

## PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Incidin™ OxyWipe S	6 x 100 wipes	3116040
Incidin™ OxyWipe S XL	6 x 50 wipes	3116700

ECOLAB HEALTHCARE  
Winnington Avenue,  
Northwich,  
Cheshire CW8 4DX  
United Kingdom



ECOLAB DEUTSCHLAND GmbH  
Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany

