INSTRUMENT
REPROCESSING



DISCOVER **NEW SYNERGIES** FOR YOUR CSSD
5 ENZYMES FOR **1 UNIQUE PRODUCT**

Aniosyme Synergy 5

UNIQUE 5-ENZYMATIC CLEANER FOR HOLISTIC MANUAL AND AUTOMATED REPROCESSING

Aniosyme Synergy 5 is the first detergent with 5-enzymes for multiple applications. Providing excellent cleaning results for manual and automated reprocessing, Aniosyme Synergy 5 allows a one product approach for the CSSD:

Aniosyme Synergy 5:

- **UNIQUE** – 5-enzyme detergent with broad-spectrum performance
- **EFFECTIVE** – Ready from 5 minutes in washer-disinfector and all qualities of water
- **SYNERGETIC** – Our first product for manual and automated reprocessing in the CSSD



ECOLAB®

www.ecolab.com

INSTRUMENT REPROCESSING

Aniosyme Synergy 5

UNIQUE 5-enzymatic cleaner for holistic manual and automated reprocessing



Aniosyme Synergy 5
5 L can

STRONG

Proven efficacy from 5 minutes onwards combined by excellent material compatibility.

ENZYMATIC

Removing all kinds of soiling by flexible dosing due to its specially developed 5 enzymatic formula.

EFFECTIVE

Lowest concentration efficacy demonstrated from 0.1% to 0.5% provides excellent cleaning performance.

PRODUCT CHARACTERISTICS

• Density	20°C	1.10 g / ml
• pH value* (concentrate)	20°C	7.89
• pH area* (solution, 1 – 5 ml/l)	20°C	6.89 - 7.29
• Conductivity (solution, 1 – 5 ml/l in DI water)	20°C	0.08 - 0.33 mS / cm

* pH value measurements under operating conditions always by means of a pH meter with temperature compensation.

PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Aniosyme Synergy 5	2 x 5 L	3094280
Aniosyme Synergy 5	1 x 10 L	3094300

APPLICATION AREA

5 enzymes cleaning product for instrumentation. Cleaning of surgical and medical equipment, endoscopes and medical devices, in soaking, ultrasonic tanks, automatic washing machine and washing tunnel.

INSTRUCTIONS FOR USE

Manual Use

- Concentrated solution for professional use
- Use at a dilution from 0.5% to 2%
- Effective from 5 minutes
- Adapt the concentration, the contact time and the temperature of the dilution water according to the degree of soiling. Brush if necessary. Rinse thoroughly
- Maximum temperature: 35°C

Automated Use

- Concentrated solution for professional use
- Use at a dilution from 0.1% to 0.5%
- Effective from 5 minutes
- Adapt the concentration, the contact time and the temperature of the dilution water according to the degree of soiling. Rinse thoroughly
- Maximum temperature: 60°C

INGREDIENTS

Non-ionic and cationic surfactants, sequestering agent, stabilising agent, enzymatic complex excipients.

STORAGE

Dangerous – respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products).
Storage: between +5°C and +25°C.

CERTIFICATES

On demand we are pleased to provide you with certificates.

MEDICAL DEVICE USE SAFELY

Please read label and product information.

+5°C +25°C

FOR PROFESSIONAL USE ONLY.

LABORATOIRES ANIOS
1, Rue de l'Espoir
59260 Lezennes
France

Distributed by:
ECOLAB HEALTHCARE
Winnington Avenue,
Northwich,
Cheshire CW8 4DX
United Kingdom

ECOLAB®

www.ecolab.com

© 2023 Ecolab. All rights reserved. HC-IR-A.SYG5-5006-2-UK-1023

0690