

VelouCid® Spray (VelouCid® Spray D)

Double-performing emulsion product to be applied post-milking by spraying.



- ▲ Powerful disinfection by iodine
- ▲ Advanced skin conditioning
- ▲ Optimum protection, in particular during cold weather



APPLICATION

VelouCid® Spray is a ready-to-use product. After each milking apply **VelouCid® Spray** on the whole teat by use of a trigger or electronic sprayer like Ecolab's JetStream and do not wipe it. Keep the animals standing or walking for 5 min to ensure biocidal effect.

FEATURES AND BENEFITS

VelouCid® Spray provides a fast disinfection against aerobic bacteria (gram positive and gram negative), yeasts, and enveloped viruses according to PT03 - Veterinary Hygiene (Disinfectants).

VelouCid® Spray forms a film which is rich in skin conditioning agents due to the EmulsioCare technology. It ensures optimum protection, in particular during cold weather.

EmulsioCare, a true oil in water emulsion, enriches and moisturizes the superficial layers of the skin. **EmulsioCare** gives to **VelouCid® Spray** remarkable dermatological qualities and soothing properties during application.



SPECIFICATIONS

VelouCid® Spray	
Physical form	Liquid
Active substance	Iodine (PVP)
pH at 20°C	3.8 - 5.0
Density at 20°C	0.98 - 1.02
Storage temperature	5°C to +25°C

LOGISTICS

Packaging	20 kg	58 kg	205 kg	975 kg
Product code	2276630	2276640	2276670	2276670
EAN code	4028162276633	4028162276640	4028162276671	4028162276674

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtstrasse 7
CH-8304 Wallisellen
www.ecolab.com