

Dairy Hygiene

All-round protection

0808 178 1010

When cleaning is about application not pressure

Contact us
for your **FREE**
on farm demo

S105 FOAM KIT

**Compact, innovative system
for foam cleaning on farms**

- ✓ All-in-one solution for foaming and rinsing
- ✓ Powered by water pressure - no need for electric supply
- ✓ Use with **Inciprop Extra**



Udder Hygiene

Product	Features	Pre or post	Barrier	Cleaning	Disinfection	Cosmetic	Ecolab Technology	Active substance	Contact time	Packaging
Oxy-foam® D	Ready to use foam	Pre	-	++	++	+++	Lactic acid foam boosted by active oxygen	Lactic acid	30 Secs	20Kg, 210Kg
VelouCid® Spray	Ready to use cosmetic spray	Pre + Post	-	-	+	+++	Povidone Iodine (PVP-I) Emulsio care	Iodine	2 mins	20Kg, 205Kg
VelouCid®	Ready to use cosmetic dip	Post	-	-	+	+++	Povidone Iodine (PVP-I) Emulsio care	Iodine	2 mins	20Kg, 200Kg
IoShield® Spray*	Ready to use sprayable barrier	Post	+	-	++	++	Povidone Iodine (PVP-I) Shield	Iodine	30 Secs	20kg, 205kg, 1000kg
Io-Shield®	Ready to use barrier dip	Post	++	-	++	++	Povidone Iodine (PVP-I) Shield	Iodine	30 Secs	20Kg, 210Kg
PhytoShield®	Ready to use barrier dip	Post	+++	-	+	++	Lactic acid boosted by natural plant extracts Shield	Lactic acid	2 mins	10Kg, 20Kg
Film-Utile® Spray **	Ready to use spray	Pre + Post	-	+	+	++	-	Lactic acid	30 Secs	20kg, 215kg, 1000kg
Film-Utile® Dip	Ready to use dip	Post	-	-	+	++	-	Lactic acid	30 Secs	20Kg, 210Kg
IoKlar® Multi **	Ready to use spray and dip	Pre OR Post	-	+	++	++	Povidone Iodine (PVP-I)	Iodine	30 Secs	20Kg, 210Kg, 1000kg
Ioklar® Superdip	Ready to use dip	Post	-	-	++	++	Povidone Iodine (PVP-I)	Iodine	30 Secs	20Kg, 200kg, 1000kg
Blu-Gard® N Spray**	Ready to use spray	Pre + Post	-	+	+	++	-	Lactic acid	30 Secs	20kG, 200kG, 1000kg
Blu-Gard® N Dip	Ready to use dip	Post	-	-	+	++	-	Lactic acid	30 Secs	20Kg, 200Kg

* For use in robotic machines ** For use in both robotic machines and ADF systems



PREMIUM Range DOUBLE ACTION / CLEANING & DISINFECTION

PERFORMANCE Range SIMPLE ACTION / DISINFECTION

Milking Machine Hygiene

Product	Features	Application	Composition	Dilution	Packaging
AseptoStar	Premium alkaline cleaner and disinfectant	Circulation and bulk tank cleaner and disinfectant	Chlorine Biocide	0.4- 0.7%	24Kg
AseptoCool	Alkaline Cleaner to be used in both hot and cold water cleaning	Circulation cleaning, bulk tanks and ice tanks	Chlorine	0.7 - 1%	25Kg/265Kg
AseptoActive	Powerful alkaline, chlorine-free cleaner	Circulation cleaning and bulk tanks	Sodium and Potassium Hydroxide	0.25%	25Kg
AseptoRobot Clean	Non chlorinated liquid alkaline cleaner for automated milking systems	Circulation cleaning for automated milking machines and bulk tanks	Sodium Hydroxide plus new generation water conditioning agents.	As per manufacturers guidelines	27Kg/280Kg
AseptoPowder	Non chlorinated powder for cleaning milking machines	Circulation cleaning and bulk tanks	Sodium Hydroxide	0.5%	20Kg
AseptoPowder Hot	Non chlorinated powder for cleaning milking machine at low temperatures day to day. Plus 8 weekly chlorinated sachets for a weekly hot wash booster	Circulation cleaning and bulk tanks	Sodium Hydroxide plus 8 chlorine sachets	0.25 - 0.5%	11Kg
HorolithStar	A blended premium acid cleaner for scale removal and cleaning of milking equipment and bulk tanks	Circulation cleaning, bulk tanks and ice tanks	Nitric, Sulphuric and Phosphoric acid	0.5%	24Kg
HorolithCool	A blended acid cleaner for scale removal and cleaning of milking equipment and bulk tanks	Circulation cleaning, bulk tanks and ice tanks	Nitric and Phosphoric acid	0.7 - 1%	24Kg/250Kg
HorolithActive	Powerful acid cleaner for milking equipment and milk tanks	Circulation cleaning and bulk tanks	Phosphoric acid plus Hydrogen Peroxide	0.25%	25Kg
HorolithRobot	A blended premium acid cleaner for scale removal and cleaning of milking equipment, including milking robots and bulk tanks	Circulation cleaning for automated milking machines and bulk tanks	Phosphoric plus Citric acid	As per manufacturers guidelines	24Kg/245Kg
HorolithPowder	A blended premium acidic powder cleaner for descaling and cleaning of milking equipment and bulk tanks. Free from Nitric acid and Phosphoric acid	Circulation cleaning and bulk tanks	Sodium Bisulphate plus Citric acid	0.5%	11Kg



Animal Housing Range

Product	Features	Application	Dilution	Packaging
Inciprop Extra	High strength Liquid foaming detergent	Animal housing	2 - 3%	22Kg
Inciprop Acid	Liquid acid foam/spray cleaner	Agriculture machinery, troughs, floors and walls	2 - 5%	25Kg
Inciprop V	Liquid cleaning agent	Vehicles, engine and trailer cleaner	1%	11Kg
Incimaxx Des N	Disinfectant by peracetic acid	Surface disinfection - DEFRA approved	0.5 - 1%	24Kg



Derma Range

Product	Features	Application	Disinfection	Active substance	Dilution	Packaging
Cloth Washing						
DermaPerfect®	Liquid cleaner and disinfectant for textiles	Washing machine	++	Peracetic acid	-	20kg
DermaSoft®	Liquid cleaner and disinfectant for textiles	Bucket	++	Peracetic acid	0.5%	20kg
DermaPowder®	Powder disinfectant to clean and disinfectant textiles	Bucket	++	Peracetic acid	0.5%	8kg
Udder Care						
DermaLotion®	Liquid cleaning soap rich in skin conditioners	Bucket and cloth	-	-	0.25%	10kg, 20kg
DermaWipes®	Rich in glycerine and Aloe Vera	Preimpregnated wipes for use before milking	-	-	-	700, 2x700 refill



Ecolab Udder Hygiene Technology

PVP iodine:

Applying the minimum yet achieving the maximum protection!

Lactic Acid:


Lactic acid boosted by active oxygen and the extracts of natural plants



Active Foam Technology

A powerful foam for cleaning and disinfection of the teats

- Cleansing
- Moisturising
- Adhesive
- High Performance



Shield Technology

A second skin that physically protects the teats.



Emulsiocare Technology

The oil fat emulsion mixed in water, the best performing cosmetic technology for the teats

Disinfection power

1 | PVP Iodine



2 | Lactic acid with active oxygen



3 | Chlorhexidine



Contact us for more information:

Stretton House Derby Road Stretton Burton-on-Trent Staffordshire DE13 0DW

0808 178 1010 | info@progiene-dairy.com | www.progiene-dairy.com