



SIMPLE. CLEAN. EFFICIENT. ▶▶▶



TOPCLIN Air Freshener

AIR FRESHENER & ODOUR NEUTRALIZER

POWERED BY
ECOLAB



DESCRIPTION:

A ready-to-use air freshener specifically designed to remove unpleasant and unwanted smells.

- Non-aerosol air freshener
- Leaves a fresh long-lasting fragrance after use



USE INSTRUCTIONS:

1. Always wear recommended PPE when handling the product
2. Use directly from the bottle as a ready to use solution
3. Make sure sprayer is turned to spray mode rather than stream to allow strong atomisation
4. Using the trigger bottle, direct and spray upwards into the corners of the area
5. According to the area size spray 2-3 times to leave a pleasant scent, do not over-use



SAFETY & STORAGE:

- Do not mix with any other product.
- For further information and safety instructions, please see product label and Safety Data Sheet.
- Store between 0 °C and 50 °C in original closed container only.



APPLICATION:

Suitable for use in public and guest areas such as offices, guest rooms, washrooms and public toilets where unpleasant odours are present.

PRODUCT CODE: PACK SIZE:

9098100 6x750ml



ADDITIONAL RESOURCES



- **MSDS**
- **WALLCHARTS**
- **PRODUCT GUIDES**



Worldwide Headquarter
Ecolab Inc.
1 Ecolab Place, St. Paul, MN55102
1.800.35.CLEAN

European Headquarter
Ecolab Europe GmbH
Richtstr. 7
CH-8304 Wallisellen